



Part Paddle with Stirring Shaft
Part Number 07525-05
Serial Number P12 001083 – P12 001220



Test Certificate

		USP / Ph.Eur. Tolerances		
Dimensions	Paddle	Height	19.0 mm	± 0.5 mm
		Upper edge length	74.5 mm	± 0.5 mm
		Lower edge length	42.0 mm	± 1.0 mm
		Thickness	4.0 mm	± 1.0 mm
		Radius	1.2 mm	± 0.2 mm
		Radius	41.5 mm	± 1.0 mm
	Shaft	Diameter	9.4 – 10.1 mm	--
Material	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Navi Mumbai, 2023-04-11

11-04-2023 16:38:41

Moreshwar Pulekar
SOTAX Quality Management
C-64 | First Floor | Suchi Business Park | TTC Industrial area | MIDC Turbhe | Navi Mumbai 400705 | India